

# FilmArray® Pouch Hydration


## Pre-Run

- Approximately **800µL hydration solution** should be dispersed equally in wells 2-11, filling each well completely. A small volume of bubbles may be present in a well, but should not be greater than 10% of total well volume. Well 12 should not contain any liquid.
- Well 1 should fill completely with approximately **300µL sample solution**. Sample must first be added and mixed with the sample buffer. A small volume of bubbles may be present in the well, but should not be greater than 10% of total well volume.

## Post Run

- Wells 1-10 should have no liquid remaining with plungers fully plunged.
- Well 11 should be unchanged during the run. It should not be plunged and should contain the same amount of liquid as it did in its pre-run state.
- Well 12 should remain un-plunged, and may contain liquid from the array.

## Image of Pouch Pre-Run and Post Run



## TECHNICAL NOTE