

BIOMÉRIEUX


BIOFIRE® SPOTFIRE® SYSTEM



BIOFIRE SPOTFIRE

BIOMÉRIEUX

PIONEERING DIAGNOSTICS



*THIS CHANGES
EVERYTHING.*

BIOFIRE® SPOTFIRE® SYSTEM

FASTER. SMALLER. SMARTER.

What if you could spot the cause of their symptoms during the patient visit?

The innovative BIOFIRE SPOTFIRE System is built on proven BIOFIRE molecular PCR (Polymerase Chain Reaction) technology. It works with BIOFIRE SPOTFIRE Panels to rapidly provide syndromic results at the molecular level. The system can be operated by any healthcare professional anywhere – even at Point of Care.

- ✓ In-house PCR results
- ✓ Results during patient visit
- ✓ Multiple targets tested at once
- ✓ Easy to use with a streamlined workflow



EXPERIENCE THE BENEFITS OF SYNDROMIC TESTING

Syndromic testing takes the guesswork out of choosing which pathogens to test for by targeting multiple relevant pathogens simultaneously to detect infectious agents quickly and accurately.

The BIOFIRE SPOTFIRE Solution uses a syndromic PCR approach integrating sample preparation, amplification and analysis. With just one patient sample, it can help clinicians make treatment decisions rapidly and confidently, contributing to overall patient satisfaction.



One sample



One unique disposable



One system



Simultaneous detection of multiple pathogens



Syndromic result with easy-to-read report



Timely impact on patient management



DESIGNED FOR ONSITE TESTING



RAPID ANSWERS

Receive fast syndromic answers in about 15 minutes*.



INTUITIVE SYSTEM INTERFACE

Designed to be run by any healthcare professional, in any setting.



SMALL FOOTPRINT

Requires minimal benchtop space with a footprint just slightly larger than a sheet of paper.



SEAMLESS CONNECTIVITY

The system can easily connect with Electronic Health Records (EHR) and other compatible systems and software.



ADVANCED USER MANAGEMENT

Administrators can set up new operators, system competency trainings, and enable initial lot lockout, among other features.



MINIMAL MAINTENANCE

No general maintenance is needed other than periodic cleaning.

SCALABLE TO YOUR NEEDS

The BIOFIRE® SPOTFIRE® System offers high throughput and random access with the ability to run up to 4 tests at the same time.



1 MODULE
26 TESTS*



2 MODULES
52 TESTS*




3 MODULES
78 TESTS*



4 MODULES
104 TESTS*

THIS CHANGES EVERYTHING.

* Calculation based on running the BIOFIRE® SPOTFIRE® Respiratory Panels (results in ~15 minutes) in an 8-hour shift.



TECHNICAL SPECIFICATIONS



MAKE THE MOST OF YOUR BIOFIRE® SPOTFIRE® INVESTMENT

Our goal is to provide you with the best support while minimizing any disruption to your workflow.



Best in-class application and engineering services



A dedicated and skilled team committed to your success, offering comprehensive global product support, specialized training courses, connectivity assistance, and expert troubleshooting for instruments and chemistry-related issues.

Product availability varies by country. Please contact your local bioMérieux representative for details.

bioMérieux S.A.
 69280 Marcy l'Etoile • France
 Tel.: + 33 (0)4 78 87 20 00 • Fax: +33 (0)4 78 87 20 90
www.biomerieux.com

Manufactured by:
 BioFire Diagnostics, LLC
 515 Colorow Drive, Salt Lake City
 UT 84108, USA

FDA CLEARED | CLIA WAIVED | CE

CUSTOMER SUPPORT

We are committed to the rapid resolution of your technical issues.



VILINK®, our proprietary service, is a highly secure remote access solution designed to maximize uptime with timely maintenance service and automatic software and firmware updates. Please contact your local bioMérieux representative or distributor for BIOFIRE products technical assistance and support.