



# VIDAS<sup>®</sup> 3

The new generation  
empowers your lab



PIONEERING DIAGNOSTICS

## Meet any request at any time

Your first priority is to serve clinicians and their patients. The VIDAS® brand is trusted around the world as a leader in on-demand immunoassay testings.

VIDAS® 3, the new generation empowers your lab.



# VIDAS® 3

## The new generation enhances your lab

Develop your lab with **VIDAS® 3**, the new generation of VIDAS®.

### Increased performance and confidence

Automation of the analysis procedure and complete traceability create confidence and peace of mind for each result you sign.

### Greater autonomy

With access to a wide range of reagents in single-test format, you handle each request on demand, for same-day diagnostic results.

### Fast adoption and a focus on productivity

Full connectivity to internal and external networks, a short training period and minimized operator maintenance time allow your lab technicians to focus on high level tasks.

### Strength in numbers

VIDAS® gives you access to a worldwide community and a select range of services to consolidate your lab's performance and results.

## Meet any request at any time

Your first priority is to serve clinicians and their patients. The VIDAS® brand is trusted around the world as a leader in on-demand immunoassay testings.

VIDAS® 3, the new generation empowers your lab.



Height: 61 cm; 24 in / Width: 75 cm; 29.5 in / Depth: 64 cm; 25.5 in

**VIDAS® 3 is designed with your specific needs and workflow in mind to provide test results that are precise and accountable.**

The system is equipped with automation, traceability and continual access, for greater productivity.



### VIDAS® 3 allows fast adoption

- **Intuitive automation:** test loading is streamlined, ensuring standardization within your lab
- **Fully automated pipetor action** with single-use, disposable tips, to avoid contamination
- Specific **management of emergency tests**
- Any **qualified lab technician** can run routine operations



### VIDAS® 3 handles test security

- **Less human manipulation** from the sample to the result
- **Sample and reagent barcodes** are scanned upon test loading for complete traceability of all test components
- **Automated calibration, on-board dilution**
- **Integrated quality control management**
- **Bi-directional LIS connectivity**




### VIDAS® 3 focuses on productivity

- **Facilitated lab workflow** with a simple answer to complex needs
- **Reduced operator time** for maintenance, to better focus on high-level tasks
- **Remote access:** a bioMérieux specialist provides both prompt troubleshooting and data for on-site instrument maintenance via screen sharing



### VIDAS® 3 is designed for continual access

- Three separate racks each containing three configurable segments can hold a total of **27 samples, diluents and calibrators** in the temperature-controlled loading bay
- **Four independent compartments** can each run 3 tests
- VIDAS® 3 can process up to **36 tests per hour**



### VIDAS® 3 uses VIDAS® reagents

- **Over 100 diagnostic tests** in single-test format for cardiovascular and infectious diseases, cancer, fertility, thyroid and allergies
- **On-demand testing:** 1 patient, 1 test, 1 result
- **Ready-to-use reagents**
- **One order reference per test:** kits contain all components needed for one test



## The new generation is more connected

You benefit from bioMérieux's experience with VIDAS® services, ensuring your lab's continued development.

### INSTRUMENT SERVICES

- Remote support services provide solutions, immediately
- Extensive network of local subsidiaries and distributors: skilled systems engineers and application specialists available wherever you are

### CERTIFICATION SUPPORT AND QUALITY CONTROL

- VIDAS® validation package for your instrument and methodology
- Support in obtaining ISO 15189 certification
- Quality data management with integrated QC software

### ONLINE TRAINING

- Online training sessions complement VIDAS® on-site service
- Continually develop your lab's skills and expertise

### A WORLDWIDE VIDAS® COMMUNITY

- The myvidas.com platform strengthens your expertise and makes the most of your VIDAS® solution